



WWW.FYICHINA.COM.CN

FANYUAN INSTRUMENT (HF) CO., LTD.  
DOWELL SCIENCE AND TECHNOLOGY (HK) CO., LTD.  
ADD: Tianyuan Road 9#, High-tech Zone, Hefei City, China Postal Code: 230088  
TEL: +86-551-65232508  
FAX: +86-551-65232507  
Email: sales@hkdownell.com  
Website: <http://www.fyichina.com.cn>

## F600 Fiber Blend Analyser





## Model No.

**F600**

## Application:

**F600 Fibre Blend Analyser** is used to quickly and accurately analysis to be made of a two component blend of wool and synthetic fiber, using a chemical method as specified.

## Related Standards:

<b>STANDARDS</b>	<b>BS 4407—1988</b>	<b>M&amp;S P100</b>
------------------	---------------------	---------------------

Note: this tester can conform to but not limit to all the standards above, for more standards conformance, please contact us.

## Features:

- Free-standing;
- Automatically feed chemicals or distilled water to sample chamber through special peristaltic metering pump(famous brand);
- Be able to pre-heat solution via IR device
- 6-working position, to compare testing result with one another
- quick, safe reverse-flushing function to avoid being blocked by specimen while filtering
- large size LCD to show testing procedures

## Key Specifications:

- Testing range 0.1%-100%
- Specimen range 0.5-3g
- Working position 6
- Power supply 220V 50HZ 3.2KW
- Dimension 560\*480\*765mm
- Net weight 45Kg



WWW.FYICHINA.COM.CN

FANYUAN INSTRUMENT (HF) CO., LTD.  
DOWELL SCIENCE AND TECHNOLOGY (HK) CO., LTD.  
ADD: Tianyuan Road 9#, High-tech Zone, Hefei City, China Postal Code: 230088  
TEL: +86-551-65232508  
FAX: +86-551-65232507  
Email: sales@hkdownell.com  
Website: <http://www.fyichina.com.cn>

**Standard configuration:**

No.	Item	Quantity
1	Main machine	1 pcs
2	Outlet/inlet hose	1pcs each
3	Drainage hose	1pcs
4	Glass Crucible	6pcs

**Optional Accessories:**

No.	Item	Quantity
1	Electronic scale 0.1mg	
2	Blast oven	